

P-97 FLAMING/BROWN MICROPIPETTE PULLER



The **P-97** micropipette puller is the latest generation of the Flaming/Brown type puller for fabrication of micropipettes, patch pipettes and microinjection needles. While retaining many of the features of earlier models, the **P-97** offers improvements in mechanical, electronic and software design. The result is better control of the pulling process and a higher degree of reproducibility. The **P-97** combines a proven mechanical system with a sophisticated, programmable microprocessor controller. This programmable control of the pulling parameters allows the investigator to design application specific pipettes from a wide range of glass compositions and sizes.

A number of other features have been incorporated in the design of the **P-97**. One

of the most apparent is the environmental chamber which houses the area surrounding the heating filament. The environmental chamber is designed to minimize the effect of changing humidity on the reproducibility of pulled pipettes. A 25% increase in power over the previous version allows for the use of larger heating filaments, larger diameter glass and multi-barreled glass. The metal jaws that clamp the heating filament have also been redesigned to minimize heat retention. There are two modes of cooling: time and delay. The delay mode provides extended cooling for large diameter and multi-barreled glass. A spring-loaded clamping mechanism has been added for easier loading of glass. A vacuum fluorescent display has been added that allows easy viewing from any direction.

Software improvements on the **P-97** include a display of the total time that the heat is on to assist in program development and troubleshooting. Up to 100 programs can now be written and stored in memory, which makes the **P-97** suitable for multiple users. These programs can now be write-protected, adding security to prevent programs from being changed or altered inadvertently. The display shows the last date and time the program was written or edited. In addition, the air pressure is a programmable parameter.

* Patent No. 4,600,424

FEATURES

> P-97

Environmental chamber

Programmable air pressure

Memory storage for up to 100 programs

Write protection and date stamp for each program

Two symmetrical pipettes with each pull

Two cooling modes: time and delay

Preprogrammed sample programs for intracellular and patch pipettes. Special programming on request

Ramp test to establish program heat settings when a new filament or glass is introduced

Vacuum fluorescent display

Internal memory test

Constant current power supply for filament and pull solenoid

Looping pull cycle for fabrication of patch type micropipettes

Self-contained air supply with filtration system and humidity control

Consistent and reliable electrodes with tip diameters less than $0.1\mu\text{m}$

Control over the time and pressure at which the air is delivered

Optimized velocity sensing circuit for maximized sensitivity and reproducibility

Quality control, SEM photograph of a tip pulled with each puller; criterion is tip measurement less than $0.1\mu\text{m}$ and typically is $\sim 0.06\mu\text{m}$

SPECIFICATIONS**> P-97****Dimensions**

21in x 14in x 12in
53cm x 36cm x 30cm

Weight

62lbs
29kg

Electrical

115/230 Volts
50/60 Hertz power line

**REFERENCES > P-97**

These references describe the Flaming/Brown series of pullers and contain valuable information applicable to the **P-97**.

1. Brown, K.T. and Flaming, D.G. *Neurosciences Journal*: 2:813-827, 1977.
2. Flaming, D.G. and Brown, K.T., *Journal of Neuroscience Methods*: 6:91-102, 1982.
3. Brown, K.T. and Flaming, D.G., *Advanced Micropipette Techniques for Cell Physiology*. John Wiley and Sons. Great Britain, 1986.

US PRICES> **P-97**

*International prices are 5% higher.
Prices subject to change without notice.*

P-97 Flaming/Brown type micropipette puller, manual, **\$ 6,950**
and Pipette Cookbook

Each puller comes with an assortment of 4 filaments and a sample box of BF100-50-10 and BF150-86-10 glass. Sutter pre-programs the P-97 with a 2.5mm box filament unless an alternative filament is requested.

ACCESSORIES

FPS	Fire polishing spacer for P-97 and P-87 pullers	\$ 25
CTS	Ceramic tile for scoring glass (large tips 20-200 microns)	\$ 15
BX10	Pipette storage box* (holds 10) 4 3/4 x 3 5/8 x 3/4 inches	\$ 15
BX20	Pipette storage box* (holds 20) 7 x 3 5/8 x 3/4 inches	\$ 15

** Minimum order of any 2 boxes*

