

MO-10



Single-Axis Oil Hydraulic Micromanipulator

Pencil type manipulator, easily combined with three-dimensional manual manipulator.

This is a pencil type oil hydraulic manipulator designed for use with a three-dimensional manipulator. Instead of employing a large cylinder with speed reduction gear, like the MO-8, its control unit uses a general drum type. Space-efficient installation and economical price are two notable advantages.

Accessory	Dedicated electrode holder	
Movement range	Fine	10mm, Full rotation of knob 500 μ m, Minimum graduation 2 μ m
Dimensions/ Weight	Control unit	W70 X D115 X H100mm, 1.15kg
	Drive unit	W30 X D20 X H125mm, 50g

