



Controller

The controller STM 3 features a built-in three-axes-joystick for simultaneous motorized movement of the micromanipulators DC-3K or DC-3KS in multiple directions. The height adjustment of the Z-axis is activated by turning the joystick knob, simplifying the operation as the user no longer has to look up from the microscope to switch from Y-axis to Z-axis for height adjustments.

STM 3  
joystick controller for  
DC-3K / DC-3KS



MS 314 and STM 3

Specifications

Order information

STM 3 for three axes  
Part No.: 00-75-220-0804

Accessories on request

Power supply:	12 V +/-5 %
Dimensions:	200 x 150 x 85 mm
Max. output voltage in neutral:	+/- ca. 10 V
Weight:	1.300 g

Marzhäuser GmbH & Co. KG  
In der Murch 15  
35579 Wetzlar-Steindorf  
Tel.: +49 (0) 64 41/91 16-0  
Fax: +49 (0) 64 41/91 16-40  
eMail: info@marzhauser.com  
www.marzhauser.com

For more detailed information, please visit our web-site at: [www.marzhauser.com](http://www.marzhauser.com)