



NARISHIGE

UMM-3FC

Micromanipulator

This model features a coaxial coarse and fine movement mechanism on each axis, and a slim, space-saving design. A rotation mechanism is mounted, making it easy to switch between right-and left-hand use, and to exchange pipettes. UMM-3FC is attached to a $\phi 12$ mm bar and used in combination with GJ-8 (Magnetic Stand) sold separately.



- * As a shorter working distance model, UM-3FC ($\phi 10$ mm bar-mount type) is available.
- * For coarse movement only, UMM-3C / UM-3C are available. Choose an appropriate model for the intended purpose.

Specification

Movement Range	Coarse	X/Y/Z Each 24mm
	Fine	X/Y/Z Each 6mm Full rotation of knob 250 μ m
		With Rotation Mechanism
Dimensions/Weight		W143 × D92 × H158mm, 510g

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2009