



YOU-1 Manipulator Set

Simple yet sophisticated design for outstanding versatility.

This micromanipulator set combines three-dimensional coarse movement with one-dimensional fine movement and is suitable for many different applications. The three-dimensional coarse movement unit is for positional adjustment, while the fine movement unit performs the final manipulation work. This fine movement unit has an independent angle adjustment mechanism and provides fine movement in the same direction as the holder without having to tilt the three-dimensional coarse movement unit. Inclusion of a rotation mechanism makes it easy to exchange pipettes.



Specification

Movement Range	Coarse	X/Y/Z Each 20mm
	Fine	20mm Full rotation of knob 250μm
		With Rotation Mechanism
Dimensions/Weight		W200 × D230 × H250mm, 1.75kg

Specifications are subject to change without notice.

Copyright © NARISHIGE Group. All Rights Reserved. 1999-2009