

NEW!

Magnetic Microscope Stage Platforms

NEW!

ALA Scientific introduces our economical **Microscope Stage Platforms**. These platforms are sturdy and have flexible configurations. They are made of powder coated magnetic steel. Options include screw down feet for use on vibration isolation tables with threaded tops or magnetic feet for use on magnetic table tops. Adjustable side plates for manipulator mounting is optional. The center hole is a 110mm diameter circle and fits ALA's **MS-Stage Plate**. The center platform is 8" x 14" and the side platforms are 8" x 10". The threaded legs are height adjustable from 3" to 12". The center platform size and adjustable height can accommodate most upright microscopes and X-Y translators.



Microscope Mounting Plate w/ screw down feet (**MSP-PED**)



Microscope Mounting Plate w/ magnetic feet (**MSP-MAG**)



All assembly tools included